

## **Extraoral Dental Suction System**





# Removes High Volume Of Droplets & Aerosols Produced During Treatment












#### Before use

• In the progress of treatment, high volumes of Droplets and Aerosols produced by high-speed handpieces and ultrasonic scalers affect dentist's healthy and safe.



#### After use

• ADS Extraoral Dental Suction System removes the Droplets and Aerosols from the origin, killing the chance of attaching to the dentists' bodies during treatment.



#### Ultrastrong Suction Power

 Ultrastrong Iron Heart: Wisely select the motor from reliable brand, which includes the various advantage such as negative voltage, high volume, minimal noise,long-life use, as well as the suction power is the tenfold of the central negative voltage.



#### Excellent Aerodynamics Design

• Rear-mounted Exhaust Design: Raised dust avoidance, effective suction assurance, patient comfort guarantee, minimal noise, and the air-tightness of independent R&D will prevent dirty air from leaking.



#### Medical Grade Filtration System

- **HEPA Filtration System:** 3 layers HEPA filtration system, stopping virus and germs ≥0.3µm with 99.995% efficiency
- Precise Water Vapour Filtration: Ensure dry water in the inner equipment and suction outlet.
- Medical Grade UV Light Disinfectant System: Killing all the virus and germs in the filter



#### Digitized and Intelligent Technology

- Filter Life Reminder: Display the rest life and expiration replacement reminder of the filter intelligently.
- Precise Control System: 10 levels of suction power adjustment, digitized display screen.

Dual-use per equipment: 1-3 level is the mode of air purification, 4-10 level is the mode of Droplets & Aerosols remove.



### **Extraoral Dental Suction System**





Carried Support: 800 488 9708

Add: 1590 S Milliken Ave., unit A, Ontario, CA 91761